Low and high level controls

111111 Low Level Control One (1) 3.0mL vial
111113 High Level Control One (1) 3.0mL vial

Intended use

The GloCyte Low and High Level Controls are assayed hematology controls designed to monitor the performance of the GloCyte Automated Cell Counter System. Assayed parameters include total nucleated cells and erythrocytes.

Principles of the method

Established laboratory practice is to use stable controls to monitor the performance of a diagnostic test. These controls are composed of stable cellular materials that provide a means of monitoring the performance of the GloCyte System. They are sampled in the same manner as a patient sample.

Preparation

The controls are composed of stabilized human erythrocytes and leukocytes suspended in a CSF-like matrix, with 0.1% sodium azide as a preservative.

Warnings and precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the human source material used in the manufacture of this product have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C, and HIV antibodies when tested by FDA-accepted methods. Because no test method can assure that a product derived from human blood is free of infectious agents, this material should be handled as potentially infectious. When handling this material follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures. Dispose of materials according to your institution's practices in a safe and acceptable manner that is in compliance with all local, regional, national and/or international regulations.

Caution: Avoid contact with metals; sodium azide may react with lead or copper plumbing to form highly explosive metal azides; flush with copious amounts of water when pouring dilute solutions down the drain to prevent buildup.

Storage, stability, and handling

	Storage	Stability
Unopened	2-8 °C (36-46 °F)	Refer to vial label
Opened	2-8 °C (36-46 °F)	30 days

- Tightly cap the vials after use.
- Return to 2-8 °C (36-46 °F) after use.
- Always store in the original vials.
- Do not freeze.

Indications of deterioration

Visible evidence of turbidity or color in the vials, or unacceptable results, may indicate deterioration. **Do not use the product if deterioration is suspected.**

Instructions for use

1. Remove the controls from the refrigerator and allow them to come to room temperature 18-25 °C (64-77 °F) for at least 15 minutes.
 2. Mix thoroughly by inverting the control vials 10-15 times.
- Note:** To prevent cells from settling and ensure homogeneity of the controls, mix prior to each use to obtain reproducible results.
3. Treat the control as you would a patient sample. Refer to the *User's Guide*¹ for use.
 4. Immediately recap the controls after each use and return to 2-8 °C (36-46 °F).

Expected results

Verify that the lot number on the vials matches the lot numbers on the table of assay values.

Performance characteristics

Assigned values are presented as a mean and range. The mean is derived from replicate testing of the indicated lots on multiple GloCyte Systems, operated and maintained according to the manufacturer's instructions. The range is an estimate of variation between replicates.

Assay values on a new lot of controls should be confirmed before the new lot is put into routine use. The laboratory's recovered mean should be within the assay range.

Each laboratory should establish its own mean and range parameters. The laboratory range may include values outside of the assay range values.

Limitations

The performance of this product is assured only if it is properly stored and used as described in this insert. Erroneous results can occur from adverse handling and/or storage conditions, use of expired materials, or improper mixing or usage.

For sales and service

Contact your Advanced Instruments Distributor.

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© 2021 Advanced Instruments

GloCyte is a registered trademark of Advanced Instruments

¹ GloCyte Automated Cell Counter System User's Guide (p/n 111005UG).

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